FCAL Gain Calibration Using Electrons

Adesh Subedi



INDIANA UNIVERSITY College of Arts and Sciences Bloomington

Methodology

- FCAL triggered runs (1501 1525)
- Cuts to select electrons
 - track match radius < 3 cm
 - 20 cm < r < 60 cm (614 blocks)
 - Cluster timing > 0
 - Number of Cluster = 1
- Select events in red box i.e., 0.6 < E(cluster)/p(track) < 1.8 & p<2 (GeV/c)
- For every hit in a block, sum up energy of 5x5 blocks (e25) around the block. Obtain average e25/p
- Assume average e25/ p = 1
- Gain factor = 1/observed (e25/p)



Result

e25/p distribution



- Variation in average e25/p assumed due to gain mismatch
- Assume average e25/ p = 1
- Gain factor = 1/observed (e25/p)

Testing new gains

Gain calibration can be tested by looking at LED pulser data



Effect on pi0 calibration

- Use gains for 614 channels as starting values for pi0 gain calibration for entire FCAL (see plot)
- Not much of a big improvement
- Limit pi0 calibration to only 614 channels with gains from electron calibration -limited statistics?

Reconstructed Pi0 Mass (post-cuts)

