

All Data from GlueX-doc-1474

Original table from Andrei

Shower Eg(GeV)	Theta	Phi	Sum	Upstream Npe	Fig	Downstream Npe	Fig
0.06		105 gap	inner-3		120	30	50
0.06		105 gap	inner-1		90	38	40
2.5		12 center	inner-3		7500	15	17000
2.5		12 center	inner-1		5500	34	13000
Dynamic Range					63		340
					61		325

To convert to number of photons at the SiPM sensor make the following adjustments

QE(440-520)	0.19
Atten correct	1.2
Accp guide	0.5

Shower Eg(GeV)	Theta	Phi	Sum	Upstream Photons	Fig	Downstream Photons	Fig
0.06		105 gap	inner-3		316		110
0.06		105 gap	inner-1		237		88
2.5		12 center	inner-3		16447		44737
2.5		12 center	inner-1		12061		34211
Dynamic Range					52		408
					51		390